


SONOST-2000


ULTRASOUND BMD

Specification


Measurement Method	Ultrasound
Measurement Site	Calcaneus (Heel)
Measurement Time	15 sec.
Measured parameter	T-score, Z-score, BQI, BUA, SOS Ratio of T-score and Z-score
Precision & Accuracy	1%
QC Check	Daily QC phantom
Dimension	523 X 286 X 275(357) mm




SOS (Sound of Speed)




BUS (Broadband Ultrasound Attenuation)




Calcaneus BQI & T-Score




Trend History



Calcaneus Report



Trend History



FRAX