

SONOST-3000V


ULTRASOUND BMD

Specification


Measurement Method	Ultrasound
Measurement Site	Calcaneus (Heel)
Measurement Time	15 sec.
Measured parameter	T-score, Z-score, BQI, BUA, SOS Ratio of T-score and Z-score, FRAX
Precision & Accuracy	1%
QC Check	Daily QC phantom
Dimension	615 X 310 X 386(293) mm
Weight	12.6 kg
Optional	USB Keyboard, Mouse, External Printer, Monitor



SOS (Sound of Speed)




BUA (Broadband Ultrasound Attenuation)



Calcaneus BQI & T-Score




Trend History



Calcaneus Report



Trend History



FRAX