


SONOST-3000


ULTRASOUND BMD

Specification

Measurement Method	Ultrasound
Measurement Site	Calcaneus (Heel)
Measurement Time	15 sec.
Measured parameter	T-score, Z-score, BQI, BUA, SOS Ratio of T-score and Z-score, FRAX
Precision & Accuracy	1%
QC Check	Daily QC phantom
Dimension	615 X 310 X 386(293) mm
Weight	12.6 kg
Optional	USB Keyboard, Mouse, External Printer, Monitor




SOS (Sound of Speed)




BUA (Broadband Ultrasound Attenuation)




High resolution Color TFT LCD touch screen




Calcaneus BQI & T-Score



Trend History



Calcaneus Report



Calcaneus Report



FRAX